

Bion Calibrator Data Sheet

				Bion Calibrator		LOT	J602			2028.01.30				
				Biorad			1008011/1008012			2029.07.31				
제조사	검사항목	Unit	Method	Bion Cal Conc.	Biorad 1					Biorad 2				
					Mean	SD	Range			Mean	SD	Range		
Bion	ALB	g/dL	BCG	3.6	4.0	0.2	3.6	~	4.4	2.7	0.2	2.3	~	3.1
	ALP	U/L	PNPP/AMP/IFCC	248	134	10	114	~	154	525	25	475	~	575
	AMY	U/L	EPS/IFCC	193	76	4	68	~	84	478	20	438	~	518
	AST(GOT)	U/L	UV without P5P(IFCC)	111	46	4	38	~	54	255	10	235	~	275
	ALT(GPT)	U/L	UV without P5P(IFCC)	117	32	3	26	~	38	110	8	94	~	126
	BUN	mg/dL	Urease GLDH	47.9	16.7	1.5	13.7	~	19.7	47.4	3.5	40.4	~	54.4
	CA	mg/dL	Arsenazo 3	10.8	9.0	0.8	7.4	~	10.6	12.1	1.0	10.1	~	14.1
	CHO	mg/dL	COD/POD	163	235	12	211	~	259	99	6	87	~	111
	CK	U/L	UV Rate	349	128	11	107	~	149	422	20	382	~	462
	CREA	mg/dL	Jaffe	3.82	2.63	0.20	2.23	~	3.03	5.55	0.60	4.35	~	6.75
	U-CREA	mg/dL	Jaffe	4.46	2.97	0.40	2.17	~	3.77	6.28	0.60	5.08	~	7.48
	D-BIL	mg/dL	Oxidation	2.1	1.69	0.20	1.29	~	2.09	0.43	0.10	0.23	~	0.63
	GGT	U/L	Glupa-C IFCC	100	61	6	49	~	73	171	9	153	~	189
	GLU	mg/dL	Hexokinase	208	83	5	73	~	93	290	15	260	~	320
	GLU 1STEP	mg/dL	GOD/POD	204	84	5	74		94	298	15	268		328
	HDL	mg/dL	Direct	57	92	5	82	~	102	25	3	19	~	31
	IP	mg/dL	Molybdate	5.46	3.49	0.30	2.89	~	4.09	7.24	0.80	5.64	~	8.84
	LDH	U/L	LD(P>L) Tris	570	349	18	313	~	385	863	40	783	~	943
	LDL	mg/dL	Direct	66	99	7	85	~	113	54	4	46	~	62
	T-BIL	mg/dL	Oxidation	4.4	4.40	0.20	4.00	~	4.80	1.16	0.12	0.92	~	1.40
	TG	mg/dL	Enzymetic	144	168	10	148	~	188	97	5	87	~	107
	TP	g/dL	Biuret	5.7	6.9	0.35	6.2	~	7.6	4.5	0.2	4.1	~	4.9
	UA	mg/dL	효소법	5.40	4.54	0.35	3.84	~	5.24	10.00	0.7	8.6	~	11.4

JWBS	CA	mg/dL	Arsenazo 3	10.8	9.0	0.8	7.4	~	10.6	12.1	1.0	10.1	~	14.1
	CK	U/L	UV Rate	349	128	11	106	~	150	422	20	382	~	462
	IP	mg/dL	Molybdate	5.46	3.49	0.30	2.89	~	4.09	7.24	0.80	5.64	~	8.84
	LIPA	U/L	Colorimetry	121	62.5	4	54.5	~	70.5	74	7.5	59	~	89
	MG	mg/dL	Xylidyl Blue	2.9	2.2	0.3	1.6	~	2.8	4.7	0.4	3.9	~	5.5

*본 Data는 기종에 따라 오차가 발생할 수도 있음을 알려드립니다.